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Artikelnummer:
626 848

Liebe Kundin, lieber Kunde!

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
en

Massage Exercises

www.tchibo.de/instructions
(Please enter the product number in the box labelled
"Bedienungsanleitungssuche" and click on "Suchen")

Product number:
626 848

Dear Customer

Massaging your body can help ease slight aches, pains and tension in your back, arms and legs. The aim of massage is to relax and stretch the muscles and to improve blood circulation and tissue oxygenation, and by doing so to boost your sense of well-being.

This massage ball enables you to activate individual trigger points which would otherwise be hard to reach.

We hope you enjoy using your new massage ball!

Your Tchibo Team

For your safety

Read the safety warnings carefully and only use the product as described in these instructions to avoid accidental injury or damage.

Keep these instructions for future reference. If you give this product to another person, remember to include these instructions.

Intended use

The massage ball is intended for self-massage and for increasing the feeling of well-being in daily life. It is not suitable for medical or therapeutic purposes.

The massage ball is designed for private household use. It is not suitable for commercial use in fitness centres or medical/therapeutic institutions.

Important information

- Wear comfortable clothing.
- Where possible, perform the ground exercises on a hard, smooth surface that does not give way and on which the massage ball can roll easily. Thick carpets, exercise mats, etc. are generally not very suitable.
- Always inspect the massage ball before every use. If any damage is visible, discontinue use.

Consult your doctor!

- Before you start to massage yourself, consult your doctor and ask about the extent to which the exercises are suitable for you.
- Discontinue massaging immediately if you begin to feel pain or have difficulties.
- If you have special physical impairments, such as wearing a pacemaker, if you are pregnant, or if you suffer from inflammation of the joints or tendons, have water retention in your arms or legs, orthopaedic complaints or have had surgery, then performance of the massage exercises is subject to your doctor's approval.
- Stop the massage immediately and see your doctor if any of the following symptoms occur: nausea, dizziness, excessive shortness of breath or pain in the chest region. You should also stop massaging immediately if you experience pain in joints and muscles.
- Not suitable for therapeutic massages or exercise.
- In the case of vein or connective tissue weakness, or varicose or spider veins, the leg massage exercises **8+9** should only be performed with light pressure and not over the entire length of the muscle at once, but rather in smaller sections. Discuss the type and extent of the massages with your doctor.

What to consider when doing the different exercises

- During all exercises keep upright and balanced - never arch your back! When standing, keep your legs slightly bent. Tense your bottom and stomach, but be sure never to make them cramp. The exercises should make you feel relaxed.
- Carry out all exercises at an even pace.
- Keep on breathing evenly during the exercises.
- Pay attention to your body! Move your muscles over the ball in such a way that the pressure feels good. Do not use it on parts of the body where it is not comfortable. It is important to be especially careful on and around your neck!
- Only carry out your selected exercises for as long as is comfortable for you. A few minutes will suffice - do not exercise excessively.
- Many of the exercises can be carried out standing as well as lying down. Try them out to see which one feels more comfortable for you.
- Tip: Other exercise examples and information can be found on the internet if needed.

Care

- Wipe the massage ball as required with a damp cloth. Do not use any aggressive cleaning products, or sharp or abrasive objects.
- After cleaning, let the massage ball air-dry at room temperature. Do not place it on a heater and do not use a hairdryer or anything similar!
- Store the massage ball in a cool, dry place. Protect it from sunlight and pointed, sharp or rough objects and surfaces.

Material: EVA foam
Weight: approx. 26 g

1. Neck

Lean against a wall with the ball behind your neck, next to your spine.

Massage your neck muscles by moving up and down against the wall from your legs. Do not raise your shoulders when doing so.

Warning!
Do not massage directly on your spine!
Alternatively, lay on your back with the ball under your neck.

2. Shoulders and arms

With slight pressure, roll the ball over your shoulders and over the top and bottom sides of your upper arms.

Warning!
Do not massage directly on your shoulder joint!


3. Lower arm

With slight pressure, roll the ball over the top and bottom sides of your lower arms, or alternatively, lay the ball on a table and roll your lower arm over it.

4. Underneath your shoulder blades

Lay the ball under your back just underneath your shoulder blades. Find a position that is comfortable for you. Do not arch your back too much.


Move your upper body in slight circular motions. For support, stretch your arms alternately to the sides, above you, and above your head.

Warning! Do not massage directly on your spine!


5. Lower back

Lay the ball next to your spine in the lumbar region under your back. Find a position that is comfortable for you. Do not arch your back too much.

Massage your lower back by moving your upper body and hips in slight circular motions. Warning! Do not massage directly on your spine!

6. Bottom


Lay your bottom on the ball. Do not arch your back.

Massage your bottom muscles by slowly moving your hips in slight circular motions over the ball.

Alternatively, sit on the ball as illustrated.


7. Hip flexors

Position yourself with your hip flexors on the ball. Do not arch your back. Massage your hip flexors by slowly rolling your hips over the ball in slight circular motions.

8. Front thigh


Lay your front thigh on the ball. Do not arch your back.

Massage the muscles in your front thigh by slowly rolling back and forth in slight circular motions.

9. Back of thigh and calf


Lay the ball under your thigh or calf as illustrated and move your leg back and forth in slight circular motions.

10. Inner calf

With slight pressure, roll the ball over the inside of your calf.

Warning! Do not massage directly on your shinbone!

11. Foot

Place your foot on the ball while standing or sitting and slowly massage your sole in circular motions.